

Measuring Cardiac Activity: *Impedance Detection with xCELLigence RTCA Cardio System*

Introduction

iCell® Cardiomyocytes², human cardiomyocytes derived from induced pluripotent stem cells, have been optimized for rapid recovery from cryopreservation. As an extension of the validated iCell Cardiomyocytes product line, they fully recapitulate biochemical, electrophysiological, mechanical, and pathophysiological characteristics of native human cardiac myocytes. These properties combine to make iCell Cardiomyocytes² an optimal in vitro test system for interrogating cardiac biology in basic research and many areas of drug development.

The xCELLigence RTCA Cardio System (RTCA Cardio system) is a non-invasive, label-free platform that utilizes impedance changes across the cardiac monolayer to measure indirectly cardiomyocyte viability, contractility, and electrical activity. iCell Cardiomyocytes² can be cultured and maintained in an E-Plate for extended durations, thus enabling measurement of acute and sub-acute drug-induced effects. Together, iCell Cardiomyocytes² and the RTCA Cardio system offer an excellent platform for in vitro screening of compound effects on human cardiomyocyte physiology.

This Application Protocol describes how to handle iCell Cardiomyocytes² for use on the RTCA Cardio system and provides basic instructions for compound treatments, data acquisition, and analysis.

Required Equipment, Consumables, and Software

The following equipment, consumables, and software are required in addition to the materials specified in the iCell Cardiomyocytes² User's Guide.

Item	Vendor	Catalog Number
Equipment		
Multichannel Pipettor, 20 and 200 µl	Multiple Vendors	
xCELLigence RTCA Cardio System	ACEA Biosciences	
Consumables		
iCell Cardiomyocytes ² Kit (Cardiomyocytes)	Cellular Dynamics International (CDI)	CMC-100-012-000.5 (0.5 unit) CMC-100-012-001 (1 unit)
Dulbecco's Phosphate Buffered Saline without Ca ²⁺ and Mg ²⁺ (D-PBS)	Invitrogen	14190
E-Plate Cardio 96 (E-Plate)	ACEA Biosciences	06417051001 06417035001
Fibronectin	Roche Applied Science	11051407001 11080938001
Sterile Reagent Reservoirs	Multiple Vendors	


Item	Vendor	Catalog Number
Software		
RTCA Cardio Instrument Software	ACEA Biosciences	

Notes

Workflow

The cardiomyocytes are thawed and plated into an E-Plate previously coated with fibronectin. 4 hours post-plating and every 48 hours thereafter, spent medium is replaced with iCell Cardiomyocytes Maintenance Medium (Maintenance Medium). From day 4 - 8 post-plating, cells can be treated with compounds, and the cardiac activity recorded.

Note: An alternative weekend-free workflow may be acceptable. Contact CDI's Technical Support (support@cellulardynamics.com; +1 (877) 320-6688 (US toll-free) or (608) 310-5100) for more information.



Methods

Preparing the E-Plate

The E-Plate is prepared the day of plating cardiomyocytes.

1. Dilute 1 mg/ml fibronectin solution in sterile D-PBS to a final concentration of 10 µg/ml immediately before use.

Note: Reconstitute fibronectin in sterile water at 1 mg/ml according to the manufacturer's instructions. Aliquot and store at -20°C.

2. Add 50 µl/well of the 10 µg/ml fibronectin solution to the center of the wells of an E-Plate to evenly coat the bottom of the well.
3. Incubate at 37°C for at least 1 hour.

Thawing Cardiomyocytes

1. Aspirate the fibronectin solution from the E-Plate. Immediately add 50 µl/well of 37°C iCell Cardiomyocytes Plating Medium (Plating Medium) to the center of the wells.

Note: Do not allow the fibronectin-coated surface to dry.

2. Equilibrate the E-Plate in a cell culture incubator at 37°C, 5% CO₂ for 5 - 10 minutes.
3. Record a background measurement according to the RTCA Cardio Instrument Operator's Guide.

Notes

4. Thaw the cardiomyocytes according to their User's Guide to a final volume of 5 ml Plating Medium by diluting the 1 ml cell suspension from the cryovial in 1 ml of Plating Medium rinse and 3 ml of additional Plating Medium.
5. Remove a sample of cells to confirm viability using a hemocytometer (using trypan blue exclusion to identify viable cells) or an automated cell counter.
6. Calculate the final volume of Plating Medium needed to obtain a final cell plating density of 1×10^6 viable cardiomyocytes/ml using the number of viable cells/vial from the Certificate of Testing.
7. Remove the E-Plate from the RTCA Cardio Instrument and equilibrate to room temperature for 5 - 10 minutes.
8. Add 50 μ l/well of the cell suspension (50,000 cells/well) to the center of the wells using a multichannel pipettor.
9. Leave the E-Plate undisturbed in the biological safety cabinet at room temperature for 20 - 30 minutes to allow the cardiomyocytes to settle and ensure an even distribution.
10. Culture the cardiomyocytes in a cell culture incubator at 37°C, 5% CO₂ for 4 hours.

Note: Place the E-Plate in a low traffic incubator and away from the door to minimize fluctuations in temperature and air movement.

Maintaining Cardiomyocytes in the E-Plate

1. Immediately before use, equilibrate an aliquot of Maintenance Medium in a 37°C water bath.
2. Replace the Plating Medium with Maintenance Medium 4 hours post-plating. Tilt the E-Plate, remove the spent medium using a multichannel pipettor, and gently add 100 μ l/well of 37°C Maintenance Medium to the side of the well to avoid disturbing the cardiomyocyte monolayer.

Note: Do not allow the pipettor tips to touch the bottom of the well during medium removal or addition. Medium replacement may cause transient alterations to beating rhythm. Allow normal beating patterns to recover after medium replacement prior to drug application.

3. Maintain the cardiomyocytes on the E-Plate, replacing 100% of the spent medium with Maintenance Medium every 48 hours.
4. Culture the cardiomyocytes in a cell culture incubator at 37°C, 5% CO₂.
5. Perform recordings from day 4 - 8 post plating.

Data Acquisition and Analysis

The RTCA Cardio Instrument Software (RTCA software) offers a wide variety of options for data acquisition and analysis. The instructions here are meant to provide a general guidance. See the RTCA Cardio Instrument Operator's Guide for specific instructions.

The beating pattern stabilizes on day 4 post-plating cardiomyocytes into the E-Plate, when the cells can be treated with compounds and the assay can be performed to analyze the electrical activity.




Figure 1: Stably Beating iCell Cardiomyocytes² on the RTCA Cardio System
Cardiomyocytes typically beat between 30 and 65 beats per minute, with a change in amplitude greater than 0.03 cell index units and a beating rhythm irregularity less than 10% per well.

Applying Compounds

1. Immediately before use, equilibrate an aliquot of Maintenance Medium in a 37°C water bath.
2. Replace the Maintenance Medium 4 - 24 hours before recording. Tilt the E-Plate, remove the Maintenance Medium using a multichannel pipettor, and gently add 90 µl/well of Maintenance Medium to the side of the well to avoid disturbing the cardiomyocyte monolayer.

Note: *Evaporation rates can vary across the E-Plate. Changing the Maintenance Medium before compound treatment is required to ensure uniform medium volumes across the E-Plate.*

3. Culture the cardiomyocytes in a cell culture incubator at 37°C, 5% CO₂.
4. Monitor the activity of the cardiomyocytes on the E-Plate to ensure regular beating rate and stable whole-peak amplitude values are reached.
5. Prepare test compounds in Maintenance Medium at 10X the final concentration in a regular 96-well cell culture plate.

Note: *Final DMSO concentrations above 0.1% should be used with caution. Therefore, if test compounds are dissolved in DMSO, the 10X compound solutions should not exceed 1% DMSO.*

6. Equilibrate the 96-well cell culture plate containing the 10X compound solutions in a cell culture incubator at 37°C, 5% CO₂.
7. Quickly transfer 10 µl/well of the 10X compound solutions from the 96-well cell culture plate to the E-Plate. Gently mix by pipetting 3 - 5 times.

Note: *Beating rate and amplitude are temperature-dependent. The E-Plate should not be kept outside the incubator for more than 5 minutes while compounds are added.*

Data Acquisition and Analysis Using the RTCA Software

See the RTCA Cardio Instrument Software Guide for specific instructions on using the RTCA software for data acquisition and analysis.

Example Data

Beating rate, amplitude, and beating rhythm irregularity were calculated with the RTCA software. Data were normalized to the last measurement point before compound treatment and averaged across the replicate wells. Results displayed in Figures 2 and 3 were obtained with 1-minute recordings acquired 60 minutes after drug treatment.




Figure 2: iCell Cardiomyocytes² Capture Phenotypic Responses across Different Classes of Cardioactive Compounds

Modulating ion channel and GPCR activity alters the spontaneous contractile activity of iCell Cardiomyocytes². Blocking I_{Kr} , I_{Ca-L} , and I_{Na} with E4031, nifedipine, and mexiletine and stimulating the β -adrenergic pathway with isoproterenol produced the expected effects on the beat waveforms.





Figure 3: Class-specific Phenotypic Responses across Different Cardioactive Compounds Are Characterized and Quantified in iCell Cardiomyocytes² Using Selected Parameters

Blocking I_{Kr} , I_{Ca-L} , and I_{Na} with E4031 (red), mexilitine (green), and nifedipine (purple) and stimulating the β -adrenergic pathway with isoproterenol (blue) produced the expected effects on beating rate, amplitude, and rhythmicity (mean \pm SEM; $n = 3$ wells for each condition).

Summary

iCell Cardiomyocytes² provide an in vitro test system that equilibrates rapidly upon reanimation from cryopreservation to recapitulate native human cardiac myocyte physiology and function while the xCELLigence RTCA Cardio System provides a label-free technology for non-invasive monitoring of cellular functions. The methods and results presented here highlight the ease of use with which robust and relevant data can be gathered on human cardiomyocyte viability, electrical activity, and contractility. Together these tools bring 96-well based, real-time, predictive assessments of compound efficacy, potency, and toxicity on human cardiomyocytes to the drug development process.

© 2016 Cellular Dynamics International, Inc. All rights reserved.
For Life Science Research Use Only.

iCell is a registered trademark, and Cellular Dynamics and the  logo are trademarks of Cellular Dynamics International, Inc.

All other brands, product names, company names, trademarks, and service marks are the properties of their respective owners.

Version: February 2016
AP-CMC2XCL160203