

### Handling and Storage

Upon receipt, immediately transfer the cryovial to liquid nitrogen storage.

### Preparing Cell Culture Surfaces

1. Thaw laminin stock solution at room temperature.
2. Dilute laminin stock in DPBS to 10 µg/ml. Do not vortex.
3. Dispense the 10 µg/ml laminin solution into the cell culture vessel(s) according to the table below.

Culture Vessel	Volume of 10 µg/ml Laminin Solution
6-well Cell Culture Plate	1 ml
96-well Cell Culture Plate	100 µl

4. Incubate at 37°C for at least 1 hour.

### Preparing the Maintenance Medium

5. Prepare maintenance medium (see **Table 2**).
6. Filter maintenance medium using a 0.2 µm PES filter unit.
7. Store maintenance medium at 4°C for up to 2 weeks.
8. Equilibrate maintenance medium to room temperature before use.

### Thawing the Cells

1. Thaw iCell® Astrocytes cryovial in a 37°C water bath for 3 minutes. Clean with 70% ethanol.
2. Transfer the cells to a 50 ml centrifuge tube.
3. Rinse the cryovial with 1 ml of maintenance medium and add it to the centrifuge tube dropwise while swirling.
4. Slowly add 8 ml of maintenance medium to the tube dropwise while gently swirling the tube.
5. Gently mix by inverting the centrifuge tube or slowly pipetting.
6. Centrifuge at 300 x g for 5 minutes.
7. Carefully aspirate the supernatant.
8. Resuspend the cell pellet in 3 ml of maintenance medium.

### Plating the Cells


1. Remove a sample of cells to perform a cell count using a hemocytometer (using trypan blue exclusion).
2. Dilute the cell suspension with maintenance medium to obtain a desired cell plating density. The recommended cell density for most cell-based assays is 40,000 – 55,000 viable cells/cm<sup>2</sup>, but it can be assay dependent.

Culture Vessel	Surface Area (cm <sup>2</sup> )	Plating Volume	Cell Number
6-well cell culture plate	9.6	3 ml	384,000 - 528,000
96-well cell culture plate	0.32	200 µl	12,800 - 17,600

3. Aspirate the laminin solution from the pre-coated plates.
4. Dispense the cells into the cell culture vessel.
5. Culture the cells at 37°C, 5% CO<sub>2</sub>.

### Maintaining the Cells

1. Replace 50-75% of the medium every 2-3 days.
2. Culture the cells at 37°C, 5% CO<sub>2</sub>.



**Figure 1: iCell Astrocytes, 100X**  
iCell Astrocytes, 01434 at 24 hours post-plating (46,000 viable cells/cm<sup>2</sup>).

### Table 1: Required Consumables

Component	Vendor	Catalog #
Dulbecco's Phosphate Buffered Saline, No Ca <sup>2+</sup> or Mg <sup>2+</sup> (DPBS)	ThermoFisher	14190
Laminin	MilliporeSigma	L2020

### Table 2: Maintenance Medium Formulation

Component	Vendor Catalog #	Volume (ml)
DMEM/F-12, HEPES	ThermoFisher # 11330	97
Fetal Bovine Serum	GE Healthcare Life Sciences #SH30071.03	2
N-2 Supplement, 100X	ThermoFisher #17502048	1

### Contacting Technical Support

Email: [fcdi-support@fujifilm.com](mailto:fcdi-support@fujifilm.com)

Phone: 1-877-320-6688

**Conditions of Use**

The cells are for RESEARCH USE ONLY. See <https://fujifilmcdi.com/assets/tnc/standard.pdf> for USE RESTRICTIONS applicable to the cells and other terms and conditions related to the cells and their use.

**Trademarks****FUJIFILM**

iCell and MyCell are registered trademarks and the  logo is a trademark of FUJIFILM Cellular Dynamics, Inc.

All other brands, product names, company names, trademarks, and service marks are the properties of their respective owners.

**Copyright Notice**

© 2019 FUJIFILM Cellular Dynamics, Inc. All rights reserved. This document may not be reproduced, distributed, modified or publicly displayed without the prior express written permission of FUJIFILM Cellular Dynamics, Inc.

**Revision History**

Document ID: X1031

Version 1.0 August 2019

***For Research Use Only***